

# Portable Outlet UPS CPAP Battery

Rechargeable Battery with AC and USB Outlets  
and Autoswitch Technology

Item#

**CAP10082**

## Includes

- Main unit
- Wall charger/power cord
- Travel bag



## Features

- Autoswitch: when power goes out, your CPAP stays on
- Three-prong AC outlet
- Two USB ports
- LED power indication meter
- Travel bag for convenience
- Meets TSA/FAA battery guidelines
- 1-year warranty
- Rated for over 500 charge cycles
- Recharges in less than 4 hours
- May be used for other devices powered by AC or USB
- AC and USB ports may be used simultaneously

## Specs

- 2.875 lbs
- 5.375 in x 9 in x 1.5 in
- Rated Capacity: 159Wh/14.4Ah
- Rated Voltage: 14.5V
- Charge Input: 16.8V/2A
- USB Output: 5V/4.2A
- AC Output: 110v AC 60hz/100w
- AC Adapter Input: 100V-240V/50-60Hz

**NO ADDITIONAL ADAPTERS NEEDED**  
for use in the U.S.A

